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HIGH TECHNOLOGICAL LEVEL solid structure and simplicity of use make this machine suitable for high-quality average output both in the chemical - pharmaceutical and food industries.


Outer cover consists of stainless-steel panels with soundproofed inner lining and lifting-doors made of transparent polycarbonate which guarantees easy maintenance and access for format change-over and mechanical service. This machine complies with GMP guidelines as well as international safety standards. EC conformity mark.

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GUARANTEE FOR AN EFFECTIVE CUSTOMER SERVICE is our 60 years' experience in the field of powders and granulates both in the pharmaceutical and food industries.

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SET PRESS microprocessor-run electronic system detects the compression strain on each punch and the average strain of each drum turn, acting directly on the dosing unit to keep the tablet weight constant and reject the single tablets outside the pre-set tolerance limits. The printer displays optimum compression values, correction and rejection limits, count of produced and rejected tablets, date, time, etc. The system can be connected to a free-standing computer for machine remote-checking.

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
TOUCH SCREEN electronic device controlling:

- handling of machine units for dosing, compression, production and filling speed
- optional fittings: mixer, hopper product level, matrix lubrication pump, language choice, etc.

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ROTARY FEEDER two-levels/three-blades with Touch Screen adjustable speed and wide filling range allowing constant dosing of the product even at high working speeds. A gravitational feeder can be installed for granulate products.

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ASYNCHRONOUS THREE-PHASE MOTOR with inverter controlled die-holder turret speed adjustment. Machine cooling fans allow continuous processing over several shifts.

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
CENTRALISED LUBRICATION SYSTEM granting utmost protection and sliding of all mechanical parts requiring it, even under high-stress conditions. For the production of effervescent tablets, a specific matrix hole lubrication kit can be installed to facilitate tablets removal.

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